

Nitrile Examination and Protective Gloves

The instruction below should be used in conjunction with detailed information on the packaging.

Product description

Non-sterile nitrile examination and protective gloves for single use. The gloves are ambidextrous, powder-free and latex-free. The textured outer surface guarantees a good grip and the beaded cuffs make it easier to put on the gloves and prevent the gloves from slipping.

Precautions and indications for use

Wash and dry your hands carefully before taking a glove from the packaging. Before usage, inspect the gloves for any defects. Put first your dominant hand into the glove and avoid touching the working surface. Then put the next glove on your other hand by using your dominant hand already wearing the glove. After both of the gloves are on your hands adjust the fit of the gloves. Ensure by stretching that the gloves fit well against fingertips and your wrists are covered.

When taking the gloves off, start by holding on the cuff of the glove on your dominant hand with a finger of the other hand. Pull the glove carefully off inside out without touching the external surface of the glove. Then pull the other glove carefully off inside out in the same way by only touching the internal surface of the glove.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in case where glove is equal to or over 400 mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on the temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections.

Material

Nitrile

External surface

Textured

Internal surface

Polymerized / Chlorinated

Colour

Blue

Package size

100 pcs by weight

Shelf life

5 years

AQL

1.5

Intended use / Application

Intended to be used in healthcare, laboratory, cleaning, food and beauty industries to protect the user and others (e.g. patient) from cross-contamination.

Storage instructions

Store at the recommended storage temperature 5-35°C in a dry place away from sunlight and fluorescent light. Recommended relative humidity in the storage is 60 ±20%. Store the gloves in a distance of at least 1 meter from heating devices, sources of fire and ozone. Do not store the gloves in direct vicinity of solvents, oils, fuels and lubricants.

Disposal

Used gloves can be contaminated with contagious or other hazardous substances and should be disposed according to local regulation.

Components / Hazardous components

The gloves may contain components which can possibly cause irritation or allergic reactions. In case of an allergic reaction, contact healthcare professionals immediately.

26102023



Nitrile Examination and Protective Gloves

MDR classification and compliance

Class I according to Annex VIII of the Regulation (EU) 2017/745.
The gloves are compliant with the following standards:
EN 455-1:2020+A1:2022, EN 455-2:2015, EN 455-3:2015,
EN 455-4:2009, EN ISO 15223-1:2021, EN ISO 21417:2019.

PPE classification and compliance


Category III Personal Protective Equipment as per Annex I of the Regulation 2016/425. The gloves are compliant with the following standards: EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018, EN ISO 374-2:2019, EN 16523-1:2015+A1:2018, EN ISO 374-4:2019, EN ISO 374-5:2016.

Declaration of conformity and instructions for use are available at <https://www.bernermedical.fi/tuote/beat-pro-nitriilikasine-sku-bpng100/>

Manufacturer

Berner Ltd
Hitsaajankatu 24, 00810 Helsinki Finland
www.bernermedical.fi

Notified Body 2777 responsible for EU Type Examination (Module B) and Module C2
On-going Conformity:
Satra Technology Europe Ltd
Bracetown Business Park, Clonee, Dublin
D15 YN2P, Ireland



Bea Pro Nitrile glove XS

Product number: BPNG1001
Sales batch: 10 boxes
Pallet: 600 boxes
GTIN Consumer package: 6414505172980
GTIN Transport package: 6414505172997



Bea Pro Nitrile glove S

Product number: BPNG1002
Sales batch: 10 boxes
Pallet: 600 boxes
GTIN Consumer package: 6414505173000
GTIN Transport package: 6414505173017



Bea Pro Nitrile glove M

Product number: BPNG1003
Sales batch: 10 boxes
Pallet: 600 boxes
GTIN Consumer package: 6414505173024
GTIN Transport package: 6414505173031



Bea Pro Nitrile glove L

Product number: BPNG1004
Sales batch: 10 boxes
Pallet: 600 boxes
GTIN Consumer package: 6414505173048
GTIN Transport package: 6414505173055



Bea Pro Nitrile glove XL

Product number: BPNG1005
Sales batch: 10 boxes
Pallet: 600 boxes
GTIN Consumer package: 6414505173062
GTIN Transport package: 6414505173079



Test acc. to EN ISO 374-5:2016

Protection against bacteria & fungi	PASS
Protection against viruses	PASS

Test acc. to EN ISO 374-2:2019 – Level 2 (ISO 2859)

Performance level: Level 3 - AQL < 0.65
Level 2 - AQL < 1.5
Level 1 - AQL < 4.0



Test acc. to EN ISO 374-1:2016+A1:2018 (Type B)

	Level	Test acc. to EN ISO 374-4:2019 Degradation %
10% Acetic Acid	4	53.5
50% Benzalkonium Chloride*	6	29.0
4% Chlorhexidine Digluconate**	6	32.9
10% Phosphoric Acid	6	14.0
40% Sodium Hydroxide (K)	6	2.6
12% Sodium Hypochlorite	6	22.7
50% Sulphuric Acid	6	21.1
5% Ethidium Bromide	6	32.9
3% Hydrogen Peroxide	6	44.0
30% Hydrogen Peroxide (P)	2	52.8
37% Formaldehyde (T)	5	20.0
50% Glutaraldehyde	6	22.9
0.1% Phenol	6	24.7

*minimum detectable permeation rate: 5 µg/cm²/min
**minimum detectable permeation rate: 7 µg/cm²/min

Permeation performance levels as per

EN ISO 374-1:2016+A1:2018

Level 1 > 10 min
Level 2 > 30 min
Level 3 > 60 min
Level 4 > 120 min
Level 5 > 240 min
Level 6 > 480 min

EN ISO 374-4:2019 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.



Nitrile Examination and Protective Gloves

Test acc. ASTM D6978-05

Test Chemotherapy Drugs	Minimum breakthrough detection time (Specimen 1/2/3) [minutes]
Carmustine (BCNU) 3.3 mg/ml (3,300 ppm)	24.6 (24.6, 27.7, 26.9)
Cisplatin 1.0 mg/ml (1,000 ppm)	240
Cyclophosphamide (Cytoxan) 20.0 mg/ml (20,000 ppm)	240
Dacarbazine (DTIC) 10.0 mg/ml (10,000 ppm)	240
Doxorubicin Hydrochloride 2.0 mg/ml (2,000 ppm)	240
Etoposide (Toposar) 20.0 mg/ml (20,000 ppm)	240
Fluorouracil 50.0 mg/ml (50,000 ppm)	240
Methotrexate 25 mg/ml (25,000 ppm)	240
Mitomycin C 0.5 mg/ml (500 ppm)	240
Paclitaxel (Taxol) 6.0 mg/ml (6,000 ppm)	240
Thiotepa (THT) 10.0 mg/ml (10,000 ppm)	68.9 (76.6, 97.6, 68.9)
Vincristine 1.0 mg/ml (1,000 ppm)	240

Food contact

Gloves are marked with food contact symbol and comply with the requirements of Regulation (EU) No 10/2011, European Regulation (EC) No 1935/2004 and with Regulation (EC) No 2023/2006 on Good Manufacturing Practice. Gloves are suitable for handling any type of food and have been tested for Overall Migration Test acc. EN 1186.

For more detailed information please contact the manufacturer.



Nitrile Examination and Protective Gloves

SYMBOLS USED ON THE PACKAGING



Medical device



Suitable for food contact
(for details check the instruction
for use)



Personal Protective Equipment



Package made from paper, qualify
for recycling



Latex-Free



Consult instructions for use



For single use only



LOT number



Non-sterile



Expiry date



Do not use if package is damaged



Country of manufacture
accompanied by date of
manufacture



Keep away from sunlight




Manufacturer



Keep dry
(protect from moisture)



Article number



Store at the temperature 5-35°C



Package is treated as
municipal waste

ISO 374-1/TYP B



Protection against chemicals
- Marking of Type B chemical
resistant glove



Protection against chemotherapy
drugs acc. to ASTM D6978-05

ISO 374-5:2016



Protection against
micro-organism hazard (viruses,
bacteria, fungi) - acc. to EN 374-5