Genos mini

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Internal distractor for lengthening
of the metacarpalia / metatarsalia

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Internal callus distraction for lengthening of the metacarpalia / metatarsalia

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Introduction

Thanks to Ilizarov, callus distraction has become an established procedure for lengthening and axis correction of extremities. Using an external fixator that must be left in place throughout the healing phase represents a complicated and troublesome therapy for patients. To increase patient convenience and practicability and prevent pin infections at the same time, we have developed a special mini-fixator. This internal distractor allows the use of callus distraction for correcting infantile malformations (brachymetatarsia and brachymetacarpia) and treating posttraumatic conditions. It offers high stability during the distraction and knitting phases and allows bone lengthening of up to 23 mm. Thanks to its low profile and internal application, this mini-fixator provides an excellent alternative to callus distraction by means of an external fixator.

Indications

- Reconstruction of metacarpals and metatarsals
- Brachymetatarsia and brachymetacarpia

Contraindications for the internal distractor include:

- Cases of in-adequate bone volume for safe fixation of the distractor
- Open and infected wounds
- Patients with immunodeficiencies
- Apparent infections
- Hypersensitivity to foreign bodies
- Suspected sensitivity to the implant materials
- Circulatory problems, systemic diseases and metabolic + disorders
- Insufficient or inadequate bone tissue
- Secondary diseases such as degenerative processes that may adversely affect the healing process
- Regions exposed to inadequately high forces or excessive weight loads
- Patients unwilling or unable to follow instructions during the postoperative phase due to their mental, neurologic or physical condition
- Osteoporosis or osteomalacia or other severe structural bone damage that prevents the stable fixation of implant components
- Bone tumors located in the implant base region
- Apparent drug or alcohol abuse

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Follow-up treatment

Metacarpal application

Plaster immobilization for 4–6 weeks using forearm cast with finger inclusion. Regular x-ray checks. From the 4th week, rehabilitation can take place on the basis of a forearm plaster shell. Alternatively, a thermoplastic splint can be used.


Metatarsal application

Plaster immobilization for 4–6 weeks using short-leg resting cast or walking cast. Thereafter, depending on radiological follow-up results, mobilization in postoperative shoe or ready-made shoe equipped with rigid insole.

Distraction phase

The distraction process is started after a postoperative latency phase (regenerate pre-stage) of 3–5 days. We recommend using a daily distraction length of 0.5 mm. This is achieved by two full clockwise turns of the activator carried out at different times of the day. One turn is equivalent to a distraction length of 0.25 mm. During the distraction phase, proper progress is continuously monitored by medical checks paying particular attention to the tensile forces acting on the patient's soft tissues as well as to sufficient callus formation (to be monitored by X-rays). If necessary, the distraction speed must be correspondingly adjusted.

Please note that excessive soft-tissue tension can lead to elastic deformation of the distractor. This should be monitored by regular X-ray examinations of the distraction progress. If the distraction process is carried on in spite of distractor deformation, this can have adverse consequences such as overload fracture of the device, screw fracture or screw tear-out.

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
Set List

26-706-00-04	Internal Finger-/Foot Distraction
Implants	
25-772-10-09 – 17-09	Set includes Cross Drive screws Ø 2.0 mm; 1 Clip = 5 screws
25-773-10-09 – 17-09	Set includes Cross Drive screws Ø 2.3 mm; 1 Clip = 5 screws
Instruments	
25-402-99-07	Screwdriver handle, CD
25-484-97-07	Cross Drive screwdriver bit 2.0 mm
26-145-03-07	Depth gauge 2.0/2.3 mm
26-150-00-07	Handle for drill bits
26-153-16-07	Drill bits 1.5 x 70 x 35 mm dental
26-153-18-07	Drill bits 1.8 x 70 x 35 mm dental
25-650-03-04	Screw length measuring clip, red, 2.0 mm
26-750-25-07	Patient screwdriver, straight; 0.25 mm
51-400-01-07	Disconnection forceps
Storage	
55-901-01-04	Tray for implants, complete
55-901-02-04	Tray for instruments, complete

KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co., Ltd.

Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
italia@klsmartin.com

Nippon Martin K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
nippon@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 505 7838
malaysia@klsmartin.com

Martin Nederland/Marned B.V.

Huizen · The Netherlands
Tel. +31 35 523 45 38
nederland@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

London · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
P.O. Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com