



Biomek **NGenius**

NEXT GENERATION LIBRARY PREP SYSTEM

IT'S NOT JUST SMART. IT'S **INGENIOUS**.



You asked for an accurate, easy to use, cost-effective, NGS sample prep system with built-in error reduction and the flexibility to run an array of library preparation kits from multiple vendors—without requiring in-depth programming knowledge, and now it's here.

ACCELERATING
answers.

 **BECKMAN
COULTER**
Life Sciences



Introducing Biomek NGeniusS System

The **Genius** is in its Total Workflow Solution



1. Flexible application selection

As part of our comprehensive service agreement, the Biomek NGeniusS System provides a wide range of complimentary demonstrated methods from a variety of reagent manufacturers, including:

DNA Sequencing

- Whole genome
- Targeted re-sequencing (Amplicon & Hybridization Capture)

RNA Sequencing

- mRNA
- Total RNA
- Targeted RNA (Amplicon & Hybridization Capture)



2. Application control


Control user access to applications with a simple authorization process (if required). Ensure that only authorized users have access to non-lab validated applications.



3. Work-from-anywhere batch set up*

Using the **Biomek NGeniusS Portal software**, set up your batch run when and where it's most convenient for you.

Whether working from home, office, or benchtop, you can program batch sizes from 4-24 samples.



4. Reagent preparation made easy

Once you select the application and enter the number of samples, the Biomek NGeniusS software creates a **Work Aid**, which includes:

- Reagent list needed for the run
- Volume requirements for each reagent
- Instructions on how to prepare and load reagents

The Work Aid saves time and effort while ensuring reproducibility during reagent preparation.

2020-10-01-04:00:01

Gene Index

10 Samples

GENERATED: 10/1/2020 11:00AM

SECTIONS:

It is released at the end of the 10-24 long incubation.

1. Hybridization, Capture, Amplification, and Purification

REAGENT PREPARATION

REAGENT	VOLUME	CONCENTRATION	AMOUNT
Hybridization Buffer	1000 µL	1000 µL	1000 µL
Hybridization Buffer 2	1000 µL	1000 µL	1000 µL
Hybridization Buffer 3	1000 µL	1000 µL	1000 µL
Hybridization Buffer 4	1000 µL	1000 µL	1000 µL
Hybridization Buffer 5	1000 µL	1000 µL	1000 µL
Hybridization Buffer 6	1000 µL	1000 µL	1000 µL
Hybridization Buffer 7	1000 µL	1000 µL	1000 µL
Hybridization Buffer 8	1000 µL	1000 µL	1000 µL
Hybridization Buffer 9	1000 µL	1000 µL	1000 µL
Hybridization Buffer 10	1000 µL	1000 µL	1000 µL
Hybridization Buffer 11	1000 µL	1000 µL	1000 µL
Hybridization Buffer 12	1000 µL	1000 µL	1000 µL
Hybridization Buffer 13	1000 µL	1000 µL	1000 µL
Hybridization Buffer 14	1000 µL	1000 µL	1000 µL
Hybridization Buffer 15	1000 µL	1000 µL	1000 µL
Hybridization Buffer 16	1000 µL	1000 µL	1000 µL
Hybridization Buffer 17	1000 µL	1000 µL	1000 µL
Hybridization Buffer 18	1000 µL	1000 µL	1000 µL
Hybridization Buffer 19	1000 µL	1000 µL	1000 µL
Hybridization Buffer 20	1000 µL	1000 µL	1000 µL
Hybridization Buffer 21	1000 µL	1000 µL	1000 µL
Hybridization Buffer 22	1000 µL	1000 µL	1000 µL
Hybridization Buffer 23	1000 µL	1000 µL	1000 µL
Hybridization Buffer 24	1000 µL	1000 µL	1000 µL
Hybridization Buffer 25	1000 µL	1000 µL	1000 µL
Hybridization Buffer 26	1000 µL	1000 µL	1000 µL
Hybridization Buffer 27	1000 µL	1000 µL	1000 µL
Hybridization Buffer 28	1000 µL	1000 µL	1000 µL
Hybridization Buffer 29	1000 µL	1000 µL	1000 µL
Hybridization Buffer 30	1000 µL	1000 µL	1000 µL
Hybridization Buffer 31	1000 µL	1000 µL	1000 µL
Hybridization Buffer 32	1000 µL	1000 µL	1000 µL
Hybridization Buffer 33	1000 µL	1000 µL	1000 µL
Hybridization Buffer 34	1000 µL	1000 µL	1000 µL
Hybridization Buffer 35	1000 µL	1000 µL	1000 µL
Hybridization Buffer 36	1000 µL	1000 µL	1000 µL
Hybridization Buffer 37	1000 µL	1000 µL	1000 µL
Hybridization Buffer 38	1000 µL	1000 µL	1000 µL
Hybridization Buffer 39	1000 µL	1000 µL	1000 µL
Hybridization Buffer 40	1000 µL	1000 µL	1000 µL
Hybridization Buffer 41	1000 µL	1000 µL	1000 µL
Hybridization Buffer 42	1000 µL	1000 µL	1000 µL
Hybridization Buffer 43	1000 µL	1000 µL	1000 µL
Hybridization Buffer 44	1000 µL	1000 µL	1000 µL
Hybridization Buffer 45	1000 µL	1000 µL	1000 µL
Hybridization Buffer 46	1000 µL	1000 µL	1000 µL
Hybridization Buffer 47	1000 µL	1000 µL	1000 µL
Hybridization Buffer 48	1000 µL	1000 µL	1000 µL
Hybridization Buffer 49	1000 µL	1000 µL	1000 µL
Hybridization Buffer 50	1000 µL	1000 µL	1000 µL
Hybridization Buffer 51	1000 µL	1000 µL	1000 µL
Hybridization Buffer 52	1000 µL	1000 µL	1000 µL
Hybridization Buffer 53	1000 µL	1000 µL	1000 µL
Hybridization Buffer 54	1000 µL	1000 µL	1000 µL
Hybridization Buffer 55	1000 µL	1000 µL	1000 µL
Hybridization Buffer 56	1000 µL	1000 µL	1000 µL
Hybridization Buffer 57	1000 µL	1000 µL	1000 µL
Hybridization Buffer 58	1000 µL	1000 µL	1000 µL
Hybridization Buffer 59	1000 µL	1000 µL	1000 µL
Hybridization Buffer 60	1000 µL	1000 µL	1000 µL
Hybridization Buffer 61	1000 µL	1000 µL	1000 µL
Hybridization Buffer 62	1000 µL	1000 µL	1000 µL
Hybridization Buffer 63	1000 µL	1000 µL	1000 µL
Hybridization Buffer 64	1000 µL	1000 µL	1000 µL
Hybridization Buffer 65	1000 µL	1000 µL	1000 µL
Hybridization Buffer 66	1000 µL	1000 µL	1000 µL
Hybridization Buffer 67	1000 µL	1000 µL	1000 µL
Hybridization Buffer 68	1000 µL	1000 µL	1000 µL
Hybridization Buffer 69	1000 µL	1000 µL	1000 µL
Hybridization Buffer 70	1000 µL	1000 µL	1000 µL
Hybridization Buffer 71	1000 µL	1000 µL	1000 µL
Hybridization Buffer 72	1000 µL	1000 µL	1000 µL
Hybridization Buffer 73	1000 µL	1000 µL	1000 µL
Hybridization Buffer 74	1000 µL	1000 µL	1000 µL
Hybridization Buffer 75	1000 µL	1000 µL	1000 µL
Hybridization Buffer 76	1000 µL	1000 µL	1000 µL
Hybridization Buffer 77	1000 µL	1000 µL	1000 µL
Hybridization Buffer 78	1000 µL	1000 µL	1000 µL
Hybridization Buffer 79	1000 µL	1000 µL	1000 µL
Hybridization Buffer 80	1000 µL	1000 µL	1000 µL
Hybridization Buffer 81	1000 µL	1000 µL	1000 µL
Hybridization Buffer 82	1000 µL	1000 µL	1000 µL
Hybridization Buffer 83	1000 µL	1000 µL	1000 µL
Hybridization Buffer 84	1000 µL	1000 µL	1000 µL
Hybridization Buffer 85	1000 µL	1000 µL	1000 µL
Hybridization Buffer 86	1000 µL	1000 µL	1000 µL
Hybridization Buffer 87	1000 µL	1000 µL	1000 µL
Hybridization Buffer 88	1000 µL	1000 µL	1000 µL
Hybridization Buffer 89	1000 µL	1000 µL	1000 µL
Hybridization Buffer 90	1000 µL	1000 µL	1000 µL
Hybridization Buffer 91	1000 µL	1000 µL	1000 µL
Hybridization Buffer 92	1000 µL	1000 µL	1000 µL
Hybridization Buffer 93	1000 µL	1000 µL	1000 µL
Hybridization Buffer 94	1000 µL	1000 µL	1000 µL
Hybridization Buffer 95	1000 µL	1000 µL	1000 µL
Hybridization Buffer 96	1000 µL	1000 µL	1000 µL
Hybridization Buffer 97	1000 µL	1000 µL	1000 µL
Hybridization Buffer 98	1000 µL	1000 µL	1000 µL
Hybridization Buffer 99	1000 µL	1000 µL	1000 µL
Hybridization Buffer 100	1000 µL	1000 µL	1000 µL



5. Error-reduced instrument setup

Sophisticated optical analytics combined with the head-up display give you real-time feedback on labware placement, virtually eliminating loading errors.

Now you can allow less-experienced operators to run complex applications with confidence.



6. Remote system monitoring*

The 360° multicolored, pattern-encoded light bar shows current instrument status.

Not within visual range? The Biomek NGenius Portal software allows you to monitor the system status from any computer.



7. Simplified sample processing

The **Multi-functional Reaction Vessel (RV)** combines three consumables into one, reducing time spent on labware inventory management:

- Sample input plate
- Thermal cycling plate
- Clean-up plate

Track your samples—with up to **100 million unique RV barcodes**—throughout the application run.



* Requires a Google Chrome or Microsoft Edge enabled computer.

Biomek NGenius System

The latest in a long line of Smart Automation Technology

Inspired by Beckman Coulter's 30-plus years of laboratory automation and Smart Automation Technology, Biomek NGenius System represents an evolutionary leap in NGS liquid-handling systems with its best-in-class performance, error-reduction technology, ease-of-use, reduced hands-on-time, and flexibility.




Error reduction

- 1 **Dynamic DeckOptix** system provides software-guided labware placement and confirmation
- 2 **Head-up Display** ensures mission critical information is in sight
- 3 **Reagent Carousel Auto ID** confirms the correct reagents are on board for the samples being processed
- 4 **Multi-channel Liquid Level Sensing** aids in the detection of insufficient reagents, alerting the operator before it becomes a processing issue
- 5 **Sample Tracking** eliminates the worry of mixing up samples when running in pre- and post-PCR lab environments


Less hands-on-time

- 6 **Reagent Carousel** reduces the need to manually transfer reagents and automates most reagent kits designed for manual use
- 7 **Temperature-controlled Reagent Storage** protects reagents from thermal degradation and allows you to run end-to-end chemistries
- 8 **Integrated Thermal Cycler** extends walkaway time
- 9 **Multi-functional Pipetting Pod**
 - **Pipetting Range** - 1 to 1000 µL
 - **Selective Tip Type** - ensures the right tip for the right volume without manual intervention
 - **Select Tip** - 1 to 8 tips for pipetting efficiency
 - **Labware Transport** - enables longer walkaway time
- 10 **UV-C Light** standard on the system



Easier to use

- 11 Software Navigation Dial** makes system programming as simple as turn and select
- 12 Biomek NGenius Portal Software*** allows setup and monitoring of batch runs from wherever it's most convenient for you
- 13 Partial Tip Box Use** eliminates manual tip transfers and/or wasting partially used tip boxes
- 14 Multi-functional Reaction Vessel** provides 3-in-1 functionality and less inventory management
- 15 Multicolor Status Light** provides 360° line of sight system status



Flexible design

- ✓ **Dynamic Sample Batch Size** accommodates batch sizes of 4 to 24 samples
- ✓ **Complimentary Application Library** is designed to evolve with your changing chemistry needs
- ✓ **Network Enabled Scalability** enables you to seamlessly set-up and run multiple instruments from your Biomek NGenius portal.
- ✓ **Continuing Application Development** ensures your investment is protected for the foreseeable future

* Requires a Google Chrome or Microsoft Edge enabled computer.




The **Genius** is in its **Simplicity**

The unique instrument design and software functionality of the **Biomek NGenius** System not only make it easy to operate, but advanced programming skills are not necessary to set up an application and complete a run.


Set up and run a batch in 6 easy steps...

On the Biomek NGenius Portal software:


Step 1 Select an application



Step 2 Enter batch name and number of samples from 4 to 24



Step 3 Make any necessary setting adjustments



Step 4 Select the stop point that fits your schedule

Configure Run

Sections

Start at section

1. Normalize Samples

#	Section	Duration (hours:minutes)
1	Normalize Samples	0:30
2	Enzymatic Through Ligation	0:30
3	Amplify Libraries	0:30
2:00 4	Clean Up Libraries	0:30
5	Quality Check	Off System
6	Normalize and Pool Libraries	Off System

Step 5 Upload your Sample Data and select 'Ready to Run'


Sample Data

UPLOAD SAMPLE DATA TEMPLATE

Well	Sample_ID	Index	Initial Concentration	260:280	260:230	Pool
A1	1234	A1	100			1
B1	1235	B1	100			1
C1	1	C1	100			1
D1	four	D1	100			1
E1	five	E1	100			1
F1	ONE	F1	100			2
G1	seven	G1	100			2
H1	8	H1	100			2
A2	9	A2	100			2
B2	ten	B2	100			2
C2	eleven	C2	100			3
D2	twelve	D2	100			3
E2	thirteen	E2	100			3
F2	fourteen	F2	100			3
G2	fifteen	G2	100			3
H2	sixteen	H2	100			4
A3	seventeen	A3	100			4
B3	eighteen	B3	100			4
C3	nineteen	C3	100			4
D3	twenty	D3	100			4

On the instrument:

Step 6 Use the dial on the instrument to select the batch and you are ready to run





Biomek
NGeniusS

Biomek NGeniusS System

Small footprint, unimaginable power and performance

Part Number	Description	Power Requirements	Weight	Dimensions (W x D X H)
C62703	Biomek NGeniusS System	10 A, 100-240 V, 50/60 Hz	135 lb (61.23 kg)	73.2 cm x 78.5 cm x 128.5 cm (with door open)

Biomek NGeniusS Consumables

The most critical point of contact

The Biomek NGeniusS System is optimized to perform as a complete system, which is why only Biomek consumables from Beckman Coulter Life Sciences are validated and approved for use.

Part Number	Description	Quantity
C62705	Biomek Reaction Vessel, 24 well	64/Box
C62706	Biomek Lid, 24 well	96/Box
C62707	Biomek Bulk Reservoir, 20 mL/section	32/Box
C59585	Case, Tips, 1025 µL, Conductive, Filtered, 480 Tips	5 Racks/Box (96 per rack)
C62712	Case, Tips, 70 µL, Conductive, Filtered, 3840 Tips	10 Racks/Box (384 per rack)
C70672	Case, Biomek 1025 µL Tip Box, Empty, 5 Racks	5 Racks/Box
C70673	Case, Biomek 70 µL Tip Box, Empty, 10 Racks	10 Racks/Box
C70665	Biomek Seal Plate, 24 well	20/Box



© 2021 Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein (including Biomek and Biomek NGeniusS) are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

All other product names and brands are properties of their respective owners.

Biomek Automated Workstations are not intended or validated for use in the diagnosis of disease or other conditions.

In development, performance characteristics have not been validated.

For Beckman Coulter's worldwide office locations and phone numbers, please visit "Contact Us" at [beckman.com](https://www.beckman.com)

20.12.476.GEN