


**Precise measuring instruments  
for viscosity and surface science**

**Product range 2023 / 2024**

Version 2023-09-01

# Viscothermostats



Precise data. Constant values.

## Viscothermostats

LAUDA Scientific viscothermostats provide accurate tempering for all glass capillary viscometers. Various products are available to support different system configurations where international norm standards require a precise temperature stability (**ASTM 0445**, **ISO 3105** and DIN 51562 from now on **DIN 53000**).

NEW

**Control head for viscothermostats Viscocool 6, Viscotemp 18, Viscotemp ET and Viscotemp 15 / 24**

### Special Features:

- || Inbuilt cooling coil to connect an external circulation chiller (e.g. MC 250)
- || Colored TFT display for simultaneous indication of actual and set temperature
- || Powerful heater (up to 2,6 kW)



**LCD S0293**  
**Viscocool 6**  
(Image similar)

## Viscocool 6, Viscotemp 18, Viscotemp ET

For R&D applications close to ambient temperature, the **Viscocool 6** is the perfect selection. With a footprint of only 206x415 mm and an integrated thermoelectric unit it provides cooling without any additional devices.

The **Viscotemp ET** is a very economic bath made from polycarbonate where a maximum of three glass capillary viscometers for manual measurement using the stop watch or up to two automatic measuring stations for iVisc or S5 (PVS) can be used. If required, when using dilution viscometers for example, up to two magnetic stirrers can be added to the thermostat.

The **Viscotemp 18** had been designed for temperatures around 100°C, e.g. for operation with silicone oils and aggressive samples. Cover plates for either five manual measuring units or for one iVisc or S5 measuring stand plus one additional glass capillary viscometer can be provided.



**LCD S0289**  
**Viscotemp ET**  
(Image similar)



Technical features	Viscocool 6 (LCD S0293)	Viscotemp 18 (LCD S0295)	Viscotemp ET (LCD S0289)
Material bath vessel	Polycarbonate	Glass	Polycarbonate
Working temperature range	15...90°C	0 <sup>1)</sup> ...105 °C	20 <sup>1)</sup> ...100 °C
Temperature stability	± 0.01 K	± 0.01 K	± 0.01 K
Max. amount measuring stands	1	1	2
Max. amount manual thermost. positions	2	5	3
VRM mounting set available	✗	✗	✓
Heater power	2.6 kW	2.6 kW	2.6 kW
Bath volume	6.5 l	18.5 l	15 l
Voltage	230 V, 50/60 Hz <sup>2)</sup>	230 V, 50/60 Hz <sup>2)</sup>	230 V, 50/60 Hz <sup>2)</sup>
Dimensions (W x D x H)	206 x 415 x 530 mm	Ø 310 x 510 mm	428 x 130 x 532 mm
Supplied cover plate	Cover plate Viscocool A1M1	Cover plate Viscotemp 18 M5	Cover plate Viscotemp ET A2

<sup>1)</sup> Can be realized with external cooling

<sup>2)</sup> Other power supply variants on request

## Viscotemp 15 / 24


The Viscotemp 15 / 24 provides a high-quality stainless steel bath and glass window and can be used in the most common temperature range from 0 up to 105 °C. The well-arranged single-chamber system with optional background lighting allows a crystal-clear view and can be easily cleaned. Typical applications for that products are the determination of the viscosity index of motor oils or the solution viscosity of plastics.

The **Viscotemp 15** can be used with up to four manual measuring stations or two automatic measuring stands, iVisc or S5 (PVS). **Viscotemp 24** offers the capacity to work with maximum seven viscometers or up to four automatic measuring stands. On both versions the VRM 4 cleaning module can be installed. The corrosion-resistant thermostating bath allows the handling with aggressive samples, e.g. polyamides dissolved in sulfuric acid and is approved for the usage of suitable silicon oils (e.g. Therm 180) as heat transfer liquid.

### Special Features:

- II Corrosion-resistant stainless steel bath with 19 to 27 liter bath volumes
- II For up to 4 automatic or 7 manual measuring stations
- II Clear design makes quick cleaning possible
- II Variopump with six pumping levels and pump flow distribution for perfect homogeneity
- II Display with standards-compliant resolution of 0.01 K, intuitive operation via cursor and soft-keys

**LCD S0299**  
**Viscotemp 24**  
(Image similar)



Technical features	Viscotemp 15 (LCD S0297)	Viscotemp 24 (LCD S0299)
Working temperature range	0 <sup>1)</sup> ...105 °C	0 <sup>1)</sup> ...105 °C
Temperature stability	± 0.01 K	± 0.01 K
Max. amount measuring stands	2	4
Max. amount of manual thermostating positions	4	7
VRM mounting set available	✓	✓
Heater power	2.6 kW	2.6 kW
Bath volume	19 l	27 l
Glass pane size (W x H)	152 x 233 mm	329 x 233 mm
Voltage	230 V, 50/60 Hz <sup>2)</sup>	230 V, 50/60 Hz <sup>2)</sup>
Dimensions (W x D x H)	532 x 233 x 552 mm	708 x 233 x 552 mm
Supplied cover plate	Cover plate Viscotemp 15 A2	Cover plate Viscotemp 24 A4

<sup>1)</sup> Can be realized with external cooling

<sup>2)</sup> Other power supply variants on request

## Viscotemp Pro 15 / 24 / 36

LAUDA Scientific Viscotemp Pro viscothermostats provide an outstanding temperature stability and homogeneity over the complete temperature range for all typical viscosity applications. As such, they are ideal for use with the fully automated PVS-S5 or iVisc viscometer and they are mandatory for the use with the VAS autosampler systems. Thanks to the double-chamber principle, a constant liquid level in the measuring room is guaranteed – even when changing the temperatures.

The special **Viscotemp Pro L models** are equipped with five layers of insulating glass. By connecting a LAUDA PRO circulation thermostat RP 250 or RP 290, they are suited for low-temperature measurements down to -40 or -60 °C.

**LCD S0278**  
**Viscotemp Pro 24**  
(Image similar)



**LCD S0276**  
**Viscotemp Pro 15**  
(Image similar)



**LCD S0302**  
**Viscotemp Pro 15**  
With ball bearing when using  
silicon oil heat transfer liquids

### Special Features:

- || Corrosion-resistant stainless steel bath with 15, 24 or 36 liter bath volume
- || Double-chamber system for constant liquid level in the measuring chamber of bath
- || Multi-glazing with optional heating avoids fogged up glass panes
- || For up to 6 automatic or 9 manual measuring stations
- || Integrated heat exchanger for counter-cooling, can be combined with LAUDA circulation thermostats
- || High control precision thanks to adaptive PID regulation
- || Optional external sensor controls the temperature at the measuring location
- || Integration of magnetic stirrers and combination with VAS systems possible

Technical features	Viscotemp Pro 15 (LCD S0276)	Viscotemp Pro 24 (LCD S0278)	Viscotemp Pro 36 (LCD S0280)	Viscotemp Pro 15 L (LCD S0282)	Viscotemp Pro 24 L (LCD S0284)
Working temperature range	0 <sup>1)</sup> ...150 °C	0 <sup>1)</sup> ...150 °C	0 <sup>1)</sup> ...150 °C	-60 <sup>3)</sup> ...100 °C	-60 <sup>3)</sup> ...100 °C
Temperature stability	± 0.01 K	± 0.01 K	± 0.01 K	± 0.01 K	± 0.01 K
Max. amount measuring stands	2	4	6	2	4
Max. amount manual positions	3	7	9	3	7
Heater power	3.5 kW	3.5 kW	3.5 kW	3.5 kW	3.5 kW
Bath volume	11...15 l	19...24 l	28...36 l	11...15 l	19...24 l
Glass pane size (W x H)	149 x 230 mm	326 x 230 mm	506 x 230 mm	149 x 230 mm	326 x 230 mm
Voltage	230 V, 50/60 Hz <sup>2)</sup>	230 V, 50/60 Hz <sup>2)</sup>	230 V, 50/60 Hz <sup>2)</sup>	230 V, 50/60 Hz <sup>2)</sup>	230 V, 50/60 Hz <sup>2)</sup>
Dimensions (W x D x H)	506 x 282 x 590 mm	740 x 282 x 590 mm	1,040 x 282 x 590 mm	506 x 282 x 590 mm	740 x 282 x 590 mm

<sup>1)</sup> Can be realized with additional LAUDA cooler    <sup>3)</sup> Can be realized with LAUDA RP 290    <sup>2)</sup> Further power supply variants on request

# Cover plates for viscothermostats

## Viscotemp 18

### **LCZ 0737** Cover plate Viscotemp 18 A1M1

For 1 measuring stand and 1 manual / thermostating position

### **LCZ 0736** Cover plate Viscotemp 18 M5

For 5 manual / thermostating positions



## Viscotemp ET

### **LCZ 041** Cover plate Viscotemp ET A2


For 2 measuring stands

### **LCZ 040** Cover plate Viscotemp ET A1M2

For 1 measuring stand and 2 manual / thermostating positions

### **LCZ 0685** Cover plate Viscotemp ET M3

For 3 manual / thermostating positions



## Viscotemp 15

### **LCZ 0730** Cover plate Viscotemp 15 A2


For 2 measuring stands

### **LCZ 0731** Cover plate Viscotemp 15 A1M2

For 1 measuring stand and 2 manual / thermostating positions

### **LCZ 0729** Cover plate Viscotemp 15 M4

For 4 manual / thermostating positions



## Viscotemp 24

### **LCZ 0733** Cover plate Viscotemp 24 A4

For 4 measuring stands

### **LCZ 0734** Cover plate Viscotemp 24 A2M4

For 2 measuring stands and 4 manual / thermostating positions

### **LCZ 0732** Cover plate Viscotemp 24 M7

For 7 manual / thermostating positions



## Viscotemp Pro 15

### **LTZ 045** Cover plate Viscotemp Pro 15 A2

For 2 measuring stands

### **LTZ 048** Cover plate Viscotemp Pro 15 A1M2

For 1 measuring stand and 2 manual / thermostating positions

### **LTZ 017** Cover plate Viscotemp Pro 15 M3

For 3 manual / thermostating positions





## Cover plates for viscothermostats



### Viscotemp Pro 24

#### **LTZ 046** Cover plate Viscotemp Pro 24 A4

For 4 measuring stands



#### **LTZ 049** Cover plate Viscotemp Pro 24 A3M3

For 3 measuring stands and 3 manual / thermostating positions



#### **LTZ 019** Cover plate Viscotemp Pro 24 M5

For 5 manual / thermostating positions



#### **LTZ 023** Cover plate Viscotemp Pro 24 M7

For 7 manual / thermostating positions



### Viscotemp Pro 36

#### **LTZ 047** Cover plate Viscotemp Pro 36 A6

For 6 measuring stands



#### **LTZ 021** Cover plate Viscotemp Pro 36 M9

For 9 manual / thermostating positions



### General


#### **LTZ 052** Adapter M for cover plates A

Conversion kit for manual measurements



#### **HDR 016** Cover M complete

Use for cover plates M



#### **HDQ 090** Cover A complete

Use for cover plates A



#### **HDQ 093** Cover A for measuring stands

Use for stand position A

## Accessories for viscothermostats



#### **LCZ 9681** Heatable window frame for Viscotemp Pro 15 L

Avoids fogged up glass panes, LiBus (optional, installing by factory)

#### **LCZ 9682** Heatable window frame for Viscotemp Pro 24 L

Avoids fogged up glass panes, LiBus (optional, installing by factory)

### Bath illumination

#### **LCZ 9738** Background lighting BL 15

LED back-light panel for Viscotemp 15 and Viscotemp Pro 15 (L)

#### **LCZ 9739** Background lighting BL 24

LED back-light panel for Viscotemp 24 and Viscotemp Pro 24 (L)

#### **LCZ 9740** Background lighting BL 36

LED back-light panel for Viscotemp Pro 36

#### **EKS 068** Extension cable LiBus

5 m, for Background lighting BL 36 (LCZ 9740)

**LCZ 9738 | LCZ 9739 | LCZ 9740**

## Accessories for viscothermostats

### **LRZ 913** RS 232/485 Interface for viscothermostats

Plug-in module upgrading viscothermostat with a RS 232/485 interface, necessary for temperature control via PC with PVS software module TEMP-DLL (LDVM2023)

### **ETP 059** External temperature probe PT100-94

For viscothermostats

### **LRZ 918** Pt 100 LiBus module

Plug-in unit for use with external temperature probe PT100-94 (ETP 059)



**LRZ 913**

## Magnetic stirrer equipment for dilution viscometry

### **LMZ 841** Single-place magnetic stirrer set

For Viscotemp ET (retrofitting possible)

### **LMVZ967-1** Two-place magnetic stirrer set

For Viscotemp Pro 15 with control unit, installing by factory

### **LMVZ968-1** Four-place magnetic stirrer set

For Viscotemp Pro 24 With control unit, installing by factory

### **UD 1190** Spirit level

For viscometer alignment



**UD 1190**

## Add. external cooling equipment

### **L001046** LAUDA Microcool MC 250 (230 V / 50 Hz)

Circulation chiller, working temp. range -10...+40°C, temperature stability 0.5 ± K, other power supply variants on request, requires silicone tubing (LZS 007)

### **LMVZ884** Cool.Fix

Heat exchanger for temperatures up to -50°C, incl. 2x metal tubing MK 100 (LZM 053)

**Requires external cooling-heating circulator (e.g. LAUDA PRO RP 290 E)**

### **L000022** LAUDA PRO RP 245 E (230 V / 50 Hz)

Compact circulation thermostat for use with Viscothermostats Viscotemp Pro L 15/24, working temp. range -45...+200°C, temperature stability 0.05 ± K, heater power max. 2.5 kW, cooling output (at 20°C) 0.8 kW, (at -30°C) 0.15 kW, other power supply variants on request

### **L002502** LAUDA PRO RP 290 E (230 V / 50 Hz)

Compact circulation thermostat for use with Viscothermostats Viscotemp Pro L 15/24, working temp. range -90...+200°C, temperature stability 0.05 ± K, heater power max. 2.5 kW, cooling output (at 20°C) 0.8 kW, (at -60°C) 0.39 kW, other power supply variants on request



**LAUDA Microcool MC 250**



**LMVZ884**  
**Cool.Fix**

## Accessories for additional external cooling equipment



**LCZ 0719**  
Cooling coil set



**LZS 007**

### **LCZ 0719 Cooling coil set**

For use with circulating chillers and viscothermostats ET 15 S and Viscotemp series, also suitable for direct use with tap water

### **Tubings**

For connecting circulation chillers with the cooling coil of viscothermostats

#### **LZS 001 Silicone tubing**

8 mm inner diameter (9 mm insulated), price per meter

#### **LZS 007 Silicone tubing**

11 mm inner diameter (9 mm insulated), price per meter

#### **LZM 053 Metal tubing MK 100**

M16X1, 100 cm, for temperatures -90...150°C,  
2 pcs. needed per each thermostat

#### **LZM 055 Metal tubing MK 200**

M16X1, 200 cm, for temperatures -90...150°C,  
2 pcs. needed per each thermostat

#### **LZS 018 Viton tubing**

12 mm inner diameter (10 mm insulated), price per meter  
for temperatures -20...150°C

#### **HKA 068 Adapter**

M16X1 – M14X1,5

#### **EZS 012 Tubing clamp**

For silicon and viton tubings (LZS 001, LZS 007, LZS 018)

## Heat transfer liquids



**LZB 114 | LZB 214 | LZB 314**

### **Silicon oil**

Temperature ranges: 0...180°C<sup>1)</sup> | **10...110°C<sup>2)</sup>**

**LZB 114 Therm 180, 5 l**

**LZB 214 Therm 180, 10 l**

**LZB 314 Therm 180, 20 l**

Temperature ranges: 90...200°C<sup>1)</sup> | **90...150°C<sup>2)</sup>**

**LZB 122 Therm 250, 5 l**

**LZB 222 Therm 250, 10 l**

**LZB 322 Therm 250, 20 l**

Temperature ranges: -20...180°C<sup>1)</sup> | **-20...100°C<sup>2)</sup>**

**LZB 116 KRYO 20, 5 l**

**LZB 216 KRYO 20, 10 l**

**LZB 316 KRYO 20, 20 l**

Temperature ranges: -50...120°C<sup>1)</sup> | **-40...50°C<sup>2)</sup>**

**LZB 121 KRYO 51, 5 l**

**LZB 221 KRYO 51, 10 l**

**LZB 321 KRYO 51, 20 l**

<sup>1)</sup> Temperature range  
technically approved

<sup>2)</sup> Temperature range recommended  
for viscosity measurement



## Heat transfer liquids

### Silicon oil

Temperature ranges: -70...30 °C<sup>1)</sup> | -60...30 °C<sup>2)</sup>

**LZB 130** KRYO 95, 5 l

**LZB 230** KRYO 95, 10 l

**LZB 330** KRYO 95, 20 l

**LZB 929** Aqua Stab water stabilizer, 100 ml

Algicide for viscothermostats, not suitable for Viscocool 6 and Viscotemp ET



**LZB 929**  
Aqua Stab



[www.lauda-scientific.de](http://www.lauda-scientific.de)

**LAUDA Scientific GmbH**  
Laudaplatz 1  
97922 Lauda-Königshofen  
Germany

Phone: +49 (0) 9343 503-340  
E-Mail: [info@lauda-scientific.de](mailto:info@lauda-scientific.de)

**LAUDA**  
  
scientific