

# How Does Microchipping Work?



**1.** Scan the pet to verify whether he or she has a microchip that has been previously implanted.



**2.** Scan the Datamars blister pack to verify that the microchip is functional, and that the ID code matches the accompanying barcode labels.



**3.** Implant the Datamars microchip, then rescan the pet to verify that the microchip can be localized.



**4.** Do not forget to complete the pet's registration in a database like PetLink!



**5.** If the pet is ever lost, hopefully his or her microchip can be scanned, and the owner's information retrieved – allowing the pet to find its way home!

 petlink.net  
creating permanent bonds between owners & pets



**PETMAXX**™