

Smaller by Design

Teleflex is committed to helping make minimally invasive surgery less invasive through smaller incisions for less trauma and smaller scars.

Introducing the MiniLap Percutaneous Surgical System with MiniGrip Handle from Teleflex:

Designed by a surgeon, the slim 2.4 mm shaft diameter instruments can be percutaneously inserted using an integrated needle tip. The grasper jaws can then be deployed to grasp and manipulate tissue. The patented deployment design allows grasper jaws to open up to 12.5 mm to grasp difficult structures, like distended gall bladders, and the proprietary steel shaft provides the strength and rigidity needed for secure tissue retraction.


MiniLap can be used for a variety of procedures such as lap cholecystectomy and lap appendectomy, and benign gynecology procedures such as diagnostic laparoscopy, LSO, and endometriosis.

Fewer Trocars

- The integrated needle tip allows direct insertion through the skin, meaning fewer trocars are used.
- In many cases, 2 or more trocars can be eliminated.

Less Trauma

- The sturdy 2.4 mm shaft allows for insertion directly through the skin, leaving a nearly unnoticeable scar.
- With MiniLap, you can improve patient satisfaction while achieving the outcomes you want for your patients.



2 times larger than the Teleflex MiniLap 3 times larger than


the Teleflex MiniLap







The Patient Nearly unnoticeable scars and potentially less pain.



MiniLap System with MiniGrip Handle			
ITEM NUMBER	PRODUCT DESCRIPTION	UNIT PER BOX	
PGAC300	MiniGrip Alligator Grasper	6	
PGCC300	MiniGrip Clutch Grasper	6	
PGBK300	MiniGrip Babcock Grasper	6	
PGBC300	MiniGrip Tong Grasper	6	

MiniLap Percutaneous Surgical System				
ITEM NUMBER	PRODUCT DESCRIPTION	UNIT PER BOX		
GBC250	Alligator Grasper	6		
CLC250	Clutch Grasper	6		
BCK250	Babcock Grasper	6		
BLC250	Bowel Grasper	6		

MiniPolar Percutaneous Electrosurgical Probes

ITEM NUMBER	PRODUCT DESCRIPTION	UNIT PER BOX
ECMC300	Curved Spatula Probe	6
ECMS300	Straight Spatula Probe	6
ECMH300	Hook Probe	6
ECMP300	Conical Probe	6

To order:

Please contact us at orders.intl@teleflex.com

94 08 17 - 00 00 01 · REV A · 10 15 005 · Teleflex, the Teleflex logo, MiniGrip, MiniLap and MiniPolar are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2015 Teleflex Incorporated. All rights reserved. All data current at time of printing. Subject to technical changes without further notice.

